

(8)ALPS処理水放出口北1km

採取日	核種濃度(Bq/L ※PuはmBq/L)								
	全β放射能	¹³⁴ Cs	¹³⁷ Cs	³ H			⁹⁰ Sr	²³⁸ Pu	²³⁹⁺²⁴⁰ Pu
				減圧蒸留法	迅速分析	電解濃縮法			
R4. 5. 19	0.01	ND	0.010	ND		0.09	0.0013	ND	ND
R4. 8. 2	0.02	ND	0.004	ND		0.11	0.0009	ND	ND
R4. 11. 8	0.02	ND	0.006	ND		0.04	0.0007	ND	ND
R5. 2. 7	0.02	ND	0.007	ND		ND	ND	ND	ND
R5. 5. 10	0.01	ND	0.009	ND		ND	0.0011	ND	ND
R5. 8. 8	0.01	ND	0.007	ND		ND	ND	ND	ND
R5. 8. 25					ND				
R5. 8. 30					ND				
R5. 9. 3	0.01	ND	0.006		ND	0.08	0.0009	ND	ND
R5. 9. 12					ND				
R5. 9. 19					ND				
R5. 9. 26					ND				
R5. 10. 8					ND				
R5. 10. 12	0.01	ND	0.020		ND	0.05	0.0006	ND	ND
R5. 10. 20					ND				
R5. 10. 24					ND				
R5. 11. 3					ND				
R5. 11. 9	0.01	ND	0.011		ND	0.28	0.0012	ND	ND
R5. 11. 14					ND				
R5. 11. 22					ND				
R5. 11. 28					ND				
R5. 12. 5	0.02	ND	0.007		ND	0.06	0.0009	ND	ND
R5. 12. 15					ND				
R5. 12. 20					ND				
R6. 1. 10					ND				
R6. 1. 18	0.02	ND	0.004		ND	0.07	0.0010	ND	ND
R6. 1. 31					ND				
R6. 2. 9	0.02	ND	0.002		ND	0.05	0.0010	ND	ND
R6. 2. 15					ND				
R6. 3. 15	0.02	ND	0.009		ND	0.34	0.0005	ND	ND
R6. 4. 12	0.02	ND	0.017		ND	0.06	0.0015	ND	ND